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Số: 22/ĐÌNH BÁCH/2018

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Nguyễn Kim Tùng

Instructions for use:

After washing the breast invert the rim of the nipple shield and place the cone so that it covers the nipple and as much of the breast as possible. A small amount of lubricating jelly around the rim can be used to improve comfort and seal.

Storage & Care

- Hand wash in warm soapy water and rinse thoroughly before and after each use.
- Always sterilise before use by immersing in baby bottle sterilising solution according to the sterilising solution manufacturer's instructions, or immersing in boiling water for at least 5 minutes, or by using an electric or microwave steam steriliser.
- Some sterilising solutions can produce harmless stains.
- Store in a dry, covered container.

Caution

- Always inspect silicone breast shields before and after each use. Throw away immediately if worn, damaged or punctured.

Directions for use steriliser box:

1) Wash all accessories before sterilising. 2) Open the sterilizer box and put 25 ml of water into the base of the box. 3) Place soothers shield down or teats base down into the water. 4) Replace the steriliser box lid 5) Place the box into the microwave and microwave on full power for 90 seconds. 6) At the end of the cycle wait 3 minutes before you remove from the microwave oven.

Warning: Keep steriliser box out of the reach of children. • Microwaves can be dangerous! • Ensure you fully understand all of the features of your microwave oven before attempting to use your steriliser box. • If using a combination oven ensure the temperature of the oven is cool before placing the steriliser inside. • Also ensure that any grill features are switched off and the combination oven is working in microwave mode only.

216WW

Ref. no. 6503

© Royal Industries (Thailand) PCL 2013

MIN. 900-2500

Pur is a registered trademark in USA and EU.

www.purbrands.com
is trademark of Royal Industries (Thailand)
PCL Co., Ltd. 135 Moo 6 Sathuak 1 The Town
Klongkum In Samutprakarn 10170 Thailand
Tel: +662 2420 1333-6 Fax: +662 2420 1333
Email: info@purbrands.com
MADE IN THAILAND

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